Neuromuskulärer Stimulator

RehaBravo

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Content

Please read the Operators Manual carefully before using the device and pay attention to important references !

Page

I. Neuromuscular Stimulator RehaBravo

Introduction - What is RehaBravo ?	30
Indication	31
Contraindication	32
Advice on Application of Electrodes	33
Technical Data	34
Safety Adjustments	34
Maintenance	35
Warranty	35

II. Display Elements of the RehaBravo

Display and Operating Elements	36
Description of the Device	37

III. Application of the RehaBravo

Settings of the Preset Programmes	38
Settings of the Individual Programmes	39-40

IV. Muscle Stimulation

Types of Muscle Fibres	41-42
Treatment Effect to the Muscular Fibres	43
The Effect of Different Frequencies	44
Frequency Selection	45
Programmes 1 - 15	46-53

V. Placement of Electrodes

Electrode-Placement for Muscle-Stimulation	54-65
Electrode-Placement for TENS-Therapy	66-70

VI. Literature	71-72
-----------------------------	-------

VII. Waste Disposal	73
----------------------------------	----

Introduction – What is RehaBravo ?

RehaBravo is a neuromuscular stimulator, which has been developed according to latest findings by MTR+ in consultation with the medical profession and patients. This device is a modern 2-channel stimulator, which enables 2 pairs of electrodes (= 4 electrodes) to work simultaneously. **RehaBravo** makes high-efficient muscular stimulation possible.

The special emphasis of **RehaBravo** is the orientation to the field of orthopaedic rehabilitation. Based on years of experience in the medical field, special value was placed on careful and gentle muscle build up as well as effective regeneration by relaxation, purge and active recovery besides the numerous programmes which the device provides for the increase of muscle power and fatigue resistance. Additionally the three freely configurable programmes enable the specialist to compile specially tailored programmes for his own individual purpose. This characteristic of the **RehaBravo** offers therefore each physician or physiotherapist, but also the experienced patient, the possibility of a focused coordinated treatment for the respective state of health. As these programmes, after storing, can be changed at any time, for example to be altered according to the progress of therapy, an almost unrestricted variation range is available with this device to the user.

It should not however in this connection remain unmentioned that RehaBravo cannot achieve miracles. As far as the state of health permits, the RehaBravo should be used in consultation with the physician or physiotherapist, as far as possible as a supplement to active physical exercise.

Indication

RehaBravo was designed to aim at the following:

Muscle Stimulation:

- Avoidance and delay of muscular spasms
- Relaxation of muscular spasms
- Avoidance and delay of muscular amyotrophy
- Recovery of the muscle
- Conserve and extend the mobility
- Extension of the capillary bed
- Speeding up the process of removal of toxic substances from metabolism
- Strengthening of the muscle
- Activation of former inactive muscular fibres
- Increase of the power of resistance of the muscle against tiredness
- Avoidance of muscular imbalances

Nerve Stimulation:

- Increase of release of bodily endorphins
- Relieve or elimination of pain, for instance in the case of:
 - Cervical Spine pain
 - Thoracic spine pain
 - Lumbar spine pain
 - Sporting accidents
 - Circulation problems
 - Neuralgia
 - Myalgia
 - Scar and phantom pain
 - Arthritis
 - Fracture pain

Contraindication

RehaBravo may under no circumstances be used in the following cases:

Patients with a pace maker as functional disturbances can occur
Patients with myocard damages (damage to the cardiac muscle) or arrhythmia (inconsistent heart beat)
When operating machinery with a potential of danger
Patients with epilepsy

Please pay unconditional attention to the following ground rules:

Do not stimulate directly through ostesynthesis material (metal in the body)
Do not stimulate transcerebral (via the head)
Do not stimulate via the Carotis-Sinus-Nerve (carotid artery)
Do not stimulate in the area of the Sinus knot (heart)
Do not stimulate on the Bulbous (eye)
Do not stimulate in the Larynx-Pharynx-Area (larynx /pharynx)
Do not stimulate in places with skin defects

Take special care with the use of the electric stimulation:

If after an operation the healing process can be interrupted through a muscle contraction
After acute bleeding or bone fractures with a tendency to bleed
In the case of the menstruation
If in the case of damage to the perception nerves the normal sensation in the skin has been lost
If at the same time the connection is made with a surgical high frequency apparatus, this could result in burn marks underneath the electrodes.
If device operation occurs close to short waves or a microwave, discrepancies of the initial values of the muscle stimulator could occur.
The application of the electrodes nearby the thorax can increase the risk of ventricle cilia.
In the case of pregnancy, muscle stimulation should only be applied after consultation with the general practitioner.
Should skin irritation occur whilst using the **RehaBravo**, please interrupt the application and inform your doctor, so that an allergic reaction can be ruled out.
Patients with implants should stimulate only after a previous consultation of an experienced doctor.

Please, keep this device stored out of the reach of children.

Advice on Application of Electrodes

Make sure, that the **RehaBravo** is switched off (no display visible on the LCD-Display).

Connect the single poled ends of the electrodes plugs with the electrodes.

The output plugs of the electrode cables will be plugged into the channels A or B of the **RehaBravo**.

Furthermore make sure, that all plugs fit completely into the electrode sockets as well as the device sockets.

Thereafter place the electrodes onto the skin. An adhesive skin gel has been applied to the electrodes, which can be re-used after removal from the skin. In the case of the skin gel loosing its adhesion, please replace the electrodes.

Electrodes which are damaged due to careless handling may not be used again. Electrodes with damaged cable insulation (e.g. visible copper wire) may not be used.

Electrodes with damaged, missing or stained aluminium coating or carbon layer, also with insufficient adhesive gel may not be used.

Electrodes may not be altered or tampered with (e.g. by scissors or other tools).

Do not pull the electrodes by force on the electrode cable. To remove the electrodes from the skin, please lift the electrode at the rim towards the middle. Do not use electrodes with ointments or creams (especially which enhance the blood flow or act as local anaesthesia).

Apply the electrodes on clean and unwounded skin only.

In the case of strong bodily hair growth, the hair should be shortened with scissors.

Neither stretch the skin nor the electrodes when applying the electrodes.

Effective current densities of more than 2 mA/cm² require an increased attention of the user with regard to the intensity adjustment. All customary electrodes which are allowed for electrical stimulation can be used, if the connections match the cables of **RehaBravo**. We recommend the use of our self-adhesive electrodes with a minimum area of 40x40 mm.

Technical Data

Channels:	Two isolated independent channels
Waveform:	Asymmetrical, rectangular, bi-phasic with zero DC current
Intensity:	0 - 90 mA
Frequency:	2 - 200 Hz
Width of impulse:	50µs - 450µs
Ramp up/down time:	0,1 seconds - 9,9 seconds
Time setting:	variable: 1 minute - 9 hours + 59 minutes
Electrical supply:	9 volt block battery
Resistor:	500 Ohm 1000 Ohm 2200 Ohm
U max:	48 V 80 V 120 V

Environmental conditions for storage & transportation: -10 to +50 degree centigrade, 0-90 % atmospheric humidity, air pressure of 700-1060 hPa

Explanation of Symbols


Attention!

Please consider
Operators Manual

0123

CE-Sign

Grade of protection IP

This device is quality controlled and fulfils the rules and regulations set by the German law for medical products (MPG) as well as the guidelines of the European Community 93/42/EEG for medical products. Therefore this device carries the CE-sign „CE 0123“.

You can find the serial number of this device in the battery compartment.

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Safety Adjustments

1. The +buttons get inactive after one minute running time of a program to avoid an unintentional increase of intensity . This lock is switched off again, only by pressing the +button three-times repeatedly. Now an increasing of intensity is possible again for another minute.
2. If the electrical circuit should be interrupted for example by losing an electrode or by a broken lead wire, the device switches off immediately, if it has been above an intensity of 10 mA.
3. If any phase comes to its end the intensity falls down to 8 mA automatically and must be regulated by the user on a pleasant measure again.

Maintenance

RehaBravo has been designed to be maintenance free. Nevertheless, please adhere to the following advice:

If needed, please clean **RehaBravo** with a mild detergent, water and a damp cloth.

Do not expose **RehaBravo** to extreme humidity and do not submerge under water.

You may clean the cables with a damp cloth. You may also from time to time treat the cables with talcum powder in order to prevent cracks and to further ongoing usage.

If over a longer period **RehaBravo** is not in use, please take the battery out of the device.

Please subject the device to technical safety control in regular intervals in accordance with the local regulations.

Warranty

MTR+ Vertriebs GmbH provides a 24 month warranty, commencing from the date of our delivery (invoice date), that in case of failure, the **RehaBravo** will be serviced or replaced free of charge. All returns must first be authorised by MTR+ Vertriebs GmbH in advance.

Excluded from this warranty are failures of the **RehaBravo** which result from excessive use, careless or improper handling as well as damage by force.

The warranty loses its lawful meaning if the producer MTR+ Vertriebs GmbH has not given its consent in writing prior to repairworks or servicing being carried out. Without this written consent the manufacturer will not return the device.

Display and Operating Elements

„ON/OFF“: By pressing the button once, the device will either be switched on or off.

„+/- buttons“: This button regulates the intensity for the channels A and B. The device also has two independently controllable output channels A and B. When pressing the respective arrow buttons the intensity will either be increased or decreased. 60 seconds after the last intensity regulation, a further increase is may not be carried out due to an electronic barrier. But if the „+“ button is pressed three times, the barrier is released and the intensity can be altered again. A reduction of the intensity is possible at any time.

Whilst configuring the programmes 16 - 18 one switches between the programmable parameters by pressing the left „+/-“, button (channel A), and through pressing the right „+/-“, button (channel B) the desired value/data for respective parameter can be set. In this way, the application time, the stimulation, the pause, rising pulse and falling pulse, as well as synchronous or asynchronous stimulation can be selected.



„PRG“: With this button the preset programmes 01 – 15 or the three free programmable programmes 16 - 18 can be selected for individual use. In addition, this button stores the current adjustment in the programming mode and ends it.

„SET“: With the SET-button one can select the different programme phases. With every pressing of this button, the phase in progress is terminated the next programme phase is started. This applies to either the running of the programmes 01 – 15 and also to the programming and the sequence of programmes 16 – 18. If one keeps the SET-button pressed for 3 seconds whilst being in programmes 16 – 18, the programming part of the programme is then entered.

„LOCK“: If you want to lock the device, please open the battery compartment on the back of the unit. Press the concealed “LOCK”-button using the end of the lead wire until you hear a double beep. On the LCD-display appears “L:” and the symbol of the watch. By pressing the “PRG”-button you set the time to zero and then the device adds up all following application times. If you should have the desire to close all other programs for the use additionally, please press the “+/-”buttons after the double beep. This way you add “P” to the display, so that you can see “L: P + watch-symbol”. After the storage of this choice with the “PRG”-button the device is locked now and adjustments cannot be changed anymore. To unlock the unit press the concealed switch once again and you will hear a single beep. The information for time in use and the average mA current used for both channels can be read. By pressing the “PRG”-button you bring the device to the Home position and the stored values are deleted.

Description of the Device

Settings of the Preset Programmes 01–15

Preparations:

Make sure that **RehaBravo** is switched off. Cleanse the area of the skin which is to be treated with clear water. Following that, make sure that all plugs are completely connected to the respective electrode sockets as well as device sockets. Place the adhesive electrodes onto the area that has been prior cleansed, according to the advices in pages 32-33 of this manual.

RehaBravo is operated by an 9 volt block battery. Put the 9 volt block battery in the battery compartment. Please ensure the correct poling in accordance with the markings. If possible use alkaline batteries of a higher standard IEC 6LR61, in order to guarantee a longer period of use. If the battery voltage is too low, in the LCD-display battery will flash. Please replace batteries immediately. Empty batteries constitute special waste and therefore have to be disposed of at an appropriate collecting point. Never dispose of batteries via domestic or household waste.

Adjustment:

1. Press once the button „**On/Off**“.

RehaBravo is now switched on. The program used last, appears on the LCD-display.

2. By pressing the programme button „**PRG**“ you can choose from the programmes 01 – 15 (pre-programmed).

3. This is followed by pressing the buttons „**+/-**“.

The selected programme now runs automatically. If the button „**+**“ is pressed repeatedly the intensity of stimulation increases. If the button „**-**“, is pressed repeatedly the intensity of stimulation is reduced. For the setting of the intensity it is vital that one feels the electric current but that pain is never incurred. **Please take into account that for safety reasons the intensity falls back at every changing of phase on the value of 8 mA and must be regulated by you on the desired standard again.** Each phase end is indicated by 5 peeptones.

4. If you wish to interrupt the set programme, then press the button „**PRG**“ once. The programme will run again when pressing the button „**+**“. In this case, the intensity must be adjusted newly, too.

5. **RehaBravo** is switched off when pressing the button „**On/Off**“ or automatically one minute after the program is finished.

Settings of the Individual Programmes 16-18

1. Switch on **RehaBravo** with “**On/off**” button.
2. Select the desired program (16-18) with “**PRG**” button.
3. Press “**SET**” button for at least three seconds. **RehaBravo** changes to the programming mode and the left watch, which indicates phase duration, flashes to the sign that this function can be changed.

In the following please use the left “+” button always then, if you want to store your choice and to change to the next step. The right “+/-” buttons are respectively meant for the change of the chosen and flashing parameter.

4. Select with the right “+/-”buttons the desired time for the first phase.
5. By pressing the left “+”button you will be able to store your choice and change to the next step. In the display now the mode flashes, which you adjust according to your wishes by pushing the right “+/-”buttons.
If you have chosen “**CONT**” or “**BURST**” you can determine frequency and after storage and change with the left +button respectively the pulse width by pushing the right “+/-”buttons.
If you have chosen “**WORK/REST**”- rhythm you can determine the following parameters one after the other by pushing the right “+/-”buttons.
 - frequency (**Hz**)
 - pulse width (**µ S**)
 - working time (**sec**)
 - rest period (**sec** = recovery time)
 - increase time of the impulse (**ramp up**)
 - way down time of the impulse (**ramp down**)
 - synchronous (**SYN**) or alternating (**ALT**) work of the channels A and B. If you have selected “**SYN**” you can program **delay**-time between the channels A and B.

If you have chosen **“MOD”** you can change the following parameters one after the other.

- low limit of frequency (**Hz down**)
- upper limit of frequency (**Hz up**)
- low limit of pulse width (**μ S down**)
- upper limit of pulse width (**μ S up**)
- low limit of intensity (the display right below shows values from 0.5 to 1.0). This means that the intensity of current can be selected between 50% (=0.5) of the original adjusted mA-value 100% (=1.0). The upper limit of the intensity, that is 100%, is always the mA value adjusted by the user originally. By the fact that the intensity, even in a modulated program, can never increase above the 100%-value it is ensured, that there will be not more current, then the user has individually regulated himself during the whole stimulation process.
- modulation time (possible choice: 5 sec. up to 60 sec, except intensity: pre-set 1 sec.)

6. By pushing the **“SET”**button, you reach the second phase of the program. You can analogously program these to the first sequence now. This way you can make up to 5 phases of your own in each of the three configurable programs. If you should want to reduce a configured program of 5 phases to less phases, you simply put the time to zero in the unwanted program sections.

7. The programming is completed by pressing the **“PRG”** button and the program shall be stored as long as you don't want to change programming.

Abbreviations:

BURST: The burst mode delivers 2 blocks of 9 impulses each per second

CONT: In the continuous mode a permanent unchanged current takes place.

MOD: In the modulated mode the frequency and/or the pulse width and/or the intensity changes in a predefined rhythm automatically.

SEC: Second

Hz: Hertz is the measurement unit for the frequency of impulses per second.

μ S: Microsecond is the measurement unit for the length of every single impulse.

MA: Milli-amp is the measurement unit for the electrical current intensity.

SYN: Synchronous function of the channels A and B

ALT: Alternating function of the channels A and B

DELAY: Delay between channels A and B

Muscular Stimulation

Types of Muscle Fibres

It is a known fact that nerves control muscles by transmission of neurological codes. This code or message, depending on the required type of muscular fibre, will be transmitted in varying frequency bands. The muscle begins to contract when it receives electrical stimulus, whether the impulse is generated by the brain or by electrical stimulation. The physiological method of neuromuscular stimulation require impulses, that are similar in condition to the naturally occurring nerve signals. Provided that the imitation of nature is as exact as possible, electrical stimulation may, if necessary, be applied over a long periods of time, without causing side effects. Muscular system of the skeleton (horizontal striped musculus) consist of numerous long and thin fibres, muscular firbers that span between tendons, by which they are linked to the bones. The respective suitable stimulus is dependant on the type of muscular fibre, that is to be reached. The main distinctions between the different fibre type is as follows: -

Type I (red muscular system)

This type of fibre is also known as ST-fibre (slow twitch fibres – slowly contracting fibres). The necessary gain of energy for the functioning of the muscles happens aerobically, that is through oxidative metabolism. The mononeuron that stimulates this fibre, has a slow transmitting speed. These type of fibres are thin and of red colour (which refers back to the existence of myoglobin molecular). Internally these fibres contain a high number of mitochondria and oxidative enzymes. Type I muscular system is extremely resistant to tiredness, as it is responsible for all types of activities of tonic nature, slow and are connected with the support of bodily posture. These slow fibres are surrounded by a dense capillary network, which allows optimal functioning of the aerobic metabolism by prolonged activity in connection with a negligible expenditure of energy. The red muscle fibres bestow stability on the bodily posture muscles and support the joint. This muscle system is very important for all endurance sports like cycling, long distance running, swimming etc.

Type IIa (white muscle system)

This type is also known as FOG-fibres (fast twitch oxidative glycolytic fibres – fast contracting fibres with oxidative glycolytic metabolism). These are stimulated by a motoneuron of the phased type that has a higher transmitting speed than tonic motoneuron.

The fibres are white because of the lack of myoglobin and have a mixed metabolic activity. They are rich in glycogen and glycolytic enzymes and also contain mitochondria enzymes, the whole metabolism is more anaerobic than the aerobic-oxidative one. These fibres are also equipped with a capillary network, that transports the necessary oxygen for the aerobic process. Fibres of type IIa are therefore capable of fast contractions, which are marked by a significant expenditure of energy, that should also be maintained for long periods of time. The fibres are therefore relatively resistant to tiredness.

Type IIb (white muscular system)

This type is also known as FG-fibres (fast twitch glycol fibres – fast contracting fibres with glycolytic metabolism). This fibre is stimulated by a phasal motoneuron with a cellular body and a very big axon that transmits impulses into the muscle at very high velocity. These fibres are white and contain a high level of glycogen and glycolytic enzymes, that enables them to release a high level of anaerobic energy. The contraction is quite fast and produces great power. Due to lack of mitochondrion, these fibres will not be able to endure long-lasting activities and fatigue easily, especially in an untrained muscle. Fibres of type IIb play a very important role in all human activities, which require an explosive expenditure of energy, e.g. sports like sprinting, weight lifting, swimming, jumping etc.

Spread of muscular fibres

Nearly all human muscles show a mixed form of different fibre types. The aforementioned types of fibres appear in diverse percentages in the muscles. The relation between fibres of Type I and Type II can vary considerably. The muscles for posture (back and stomach muscular systems) show a higher proportion of Type I fibres, whereas muscles for movement / motion (leg and arm muscular systems) usually consist mainly of Type II fibres.

Treatment Effect to the Muscular Fibres

Enhancement of physical capabilities by training or electrical stimulation can be related back to the fact, that the human body is able to adjust to a certain extent to changing environmental conditions.

If the stress to a muscle is increased beyond a threshold value, a corresponding adaptation process will ensue. The increase of muscular power is initially due to improved intro-muscular co-ordination, which means that for one movement more fibres will be activated simultaneously than before. Only in the second step will the muscular mass be increased. Whether this occurs exclusively due to hypertrophy of the muscular fibres (meaning the enlargement of the fibre diameter) or if it is the result of hyperplasia (meaning the increase of muscular fibres) is scientifically not proven. It is probable that both phenomena concurrently are responsible for the muscular growth.

The aforementioned training effects to the muscular system can be attained in a more passive way by electrical stimulation. One should not make the mistake of neglecting the physical training for the better of the electrical stimulation. A movement, as the aim of almost every sporting activity, consists of a finely tuned response of different muscles or groups of muscles. Flowing movements will only be achieved through an optimum of co-ordination of the flexing of the active muscles (agonists) and relaxing of their opponents (antagonists). The flowing movements are the unalterable pre-requirements for sporting success. The therefore necessary neuromuscular co-ordination processes cannot be trained with this kind of muscle stimulation. To avoid muscular imbalances, it is urgently advisable to apply strengthening programmes for both the agonists and the antagonists in equivalent amounts.

Besides the strengthening of the muscles, is the furtherance of the muscular blood circulation an additional main aim of the electrical stimulation. A higher capillary density enable the improved blood circulation. Hence an increase supply of nutritive substances, especially oxygen, is made available and the aerobic metabolism especially benefits. Electrical stimulation not only allows an improvement of the supply, but also furthers the fast disposal of metabolic waste products like lactic acid, or carbon dioxide. The acceleration of the physical regeneration in phases of more intensive physical strain of weakened musculature is of decisive importance for a fast re-production of the normal physical capacity.

The Effect of Different Frequencies

A very short burst of electrical stimulation produces only a short contraction or „single shock“ after which the muscle falls back into its natural form and length, which it possesses in the resting position. If the stimulation occurs sequentially, we realise that the contracting effects caused by overlapping of the contracting phases are additive, as the mechanical twitches last distinctly longer than the electrical stimulation. This phenomena is called „incomplete tetanus“. Neither „single shock“ nor „incomplete tetanus“ will normally be observed in an arbitrary human movement.

The phenomena of „complete tetanus“ means the state of a muscular contraction which is caused by reoccurring electrical stimulation of the motoric nerves with a sufficiently high frequency. Hereby the individual shocks are amalgamated and are no longer distinguishable. In this situation, the muscle contracts and stiffens through the tension produced by the muscle; a measurable power is exercised at the tendon endings. Nearly all muscular contractions of human muscles possess the characteristic feature of a „complete tetanus“.

Frequency Selection

5 Hz or lower Single twitches can be achieved with this frequency, after which the muscle can relax completely for a short time. These frequencies are used in the starting phase of the strengthening programmes with the purpose of preparing the muscle for the coming working phase without tiring it. A frequency of 3 Hz lies within the frequency band for the production of bodily endorphins for the relief of pain and general relaxation.

5 – 15 Hz These frequencies are also referred to as shaking frequencies, which will not lead to the „complete tetanus“ but otherwise do not allow for complete relaxation in-between the impulses. This frequency span is selected for improvement of muscular tonus, support and stability of the joints. 10 Hz is the natural frequency of the slow oxidative muscle fibre (Type I). The electrical stimulation increases the resistance to tiredness by improving the density of the capillary system and it increases the muscles ability to cope with a reduction in oxygen. Onerous metabolic waste will be further reduced through the „milking of the veins“. This frequency band is of great importance in the regeneration phase after high levels of stress and can be used daily over a number of hours for sportive or therapeutic purpose.

15 – 30 Hz This frequency can be used for increase of stamina of the muscle. This is the natural frequency band of the oxidative-glycolytic muscle fibres (Type IIa). The frequency can be applied up to one hour daily.

30 – 60 Hz The muscle stays in a state of complete tetanus from a frequency of 30 Hz upwards. These frequencies are applied for the strengthening of the muscle and recruitment of the fast glycolytic muscle fibres (Type IIb). A treatment in these frequency bands should only be applied for relatively short periods because the tiring of the muscle will begin within a few minutes of electrical stimulation.

60 – 120 Hz These frequencies will often be selected in order that the muscle shall shed great energy (explosive energy). With this high frequency it is important to apply the stimulation only for very short periods.

Programmes

Progr.: 01		Warming up				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	6				
Mode		cont				
Frequency work	Hz	5				
Frequency rest	Hz					
Pulse duration	μ S	300				
Modulation time	secs					
Ramp up time	secs					
Ramp down time	secs					
Work time	secs					
Rest time	secs					
Alternating						
Synchronous						
Overall time	6 min					

Progr.: 02		Capillary				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	20				
Mode		Cont				
Frequency work	Hz	10				
Frequency rest	Hz					
Pulse duration	μ S	250				
Modulation time	secs					
Ramp up time	secs					
Ramp down time	secs					
Work time	secs					
Rest time	secs					
Alternating						
Synchronous						
Overall time	20 min					

Programmes

Progr.: 03		Increase of muscles 1				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	5	15	10	15	5
Mode		Cont	W/R	W/R	W/R	Cont
Frequency work	Hz	10	20	30	20	5
Frequency rest	Hz		3	3	3	
Pulse duration	μ S	250	300	300	300	250
Modulation time	secs					
Ramp up time	secs		2	2	2,5	
Ramp down time	secs		1,8	1,8	1,8	
Work time	secs		6	10	6	
Rest time	secs		10	10	10	
Alternating						
Synchronous						
Overall time	50 min					

Progr.: 04		Increase of muscles 2				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	5	15	10		
Mode		Cont	W/R	Cont		
Frequency work	Hz	5	60	2		
Frequency rest	Hz		3			
Pulse duration	μ S	300	350	250		
Modulation time	secs					
Ramp up time	secs		2			
Ramp down time	secs		1,5			
Work time	secs		7			
Rest time	secs		14			
Alternating						
Synchronous						
Overall time	30 min					

Programmes

Progr.: 05		Strengthening of muscles 1				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	5	12	5	12	6
Mode		Cont	W/R MOD	Cont	W/R MOD	Cont
Frequency work	Hz	5	50-75	5	40-75	3
Frequency rest	Hz		3		3	
Pulse duration	μ S	300	300	300	300	300
Modulation time	secs		10		10	
Ramp up time	secs		2		2	
Ramp down time	secs		1		1	
Work time	secs		10		10	
Rest time	secs		8		8	
Alternating						
Synchronous						
Overall time	40 min					

Progr.: 06		Strengthening of muscles 2				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	5	20	10		
Mode		Cont	W/R	Cont		
Frequency work	Hz	5	75	2		
Frequency rest	Hz		3			
Pulse duration	μ S	300	300	250		
Modulation time	secs					
Ramp up time	secs		1,5			
Ramp down time	secs		1			
Work time	secs		5			
Rest time	secs		12			
Alternating						
Synchronous						
Overall time	35 min					

Programmes

Progr.: 07		Severe partial peripheral paralysis				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	10				
Mode		W/R				
Frequency work	Hz	60				
Frequency rest	Hz					
Pulse duration	μ S	450				
Modulation time	secs					
Ramp up time	secs	0,6				
Ramp down time	secs	1				
Work time	secs	4				
Rest time	secs	20				
Alternating						
Synchronous						
Overall time	10 min					

Progr.: 08		Medium partial peripheral paralysis				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	15				
Mode		W/R				
Frequency work	Hz	50				
Frequency rest	Hz					
Pulse duration	μ S	300				
Modulation time	secs					
Ramp up time	secs	0,5				
Ramp down time	secs	1				
Work time	secs	5				
Rest time	secs	15				
Alternating						
Synchronous						
Overall time	15 min					

Programmes

Progr.: 09		Severe muscle atrophy				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	10				
Mode		W/R				
Frequency work	Hz	6				
Frequency rest	Hz					
Pulse duration	μ S	350				
Modulation time	secs					
Ramp up time	secs	0,8				
Ramp down time	secs	1				
Work time	secs	5				
Rest time	secs	35				
Alternating						
Synchronous						
Overall time	10 min					

Progr.: 10		Medium muscle atrophy				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	20				
Mode		W/R				
Frequency work	Hz	15				
Frequency rest	Hz					
Pulse duration	μ S	300				
Modulation time	secs					
Ramp up time	secs	0,6				
Ramp down time	secs	1				
Work time	secs	6				
Rest time	secs	30				
Alternating						
Synchronous						
Overall time	20 min					

Programmes

Progr.: 11		Scoliosis / Kyphosis				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	60				
Mode		W/R				
Frequency work	Hz	35				
Frequency rest	Hz					
Pulse duration	μ S	250				
Modulation time	secs					
Ramp up time	secs	1,5				
Ramp down time	secs	1,5				
Work time	secs	8				
Rest time	secs	20				
Alternating						
Synchronous						
Overall time	60 min					

Progr.: 12		Scoliosis 2				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	599				
Mode		WR MOD				
Frequency work	Hz	20-50				
Frequency rest	Hz					
Pulse duration	μ S	250				
Modulation time	secs	60				
Ramp up time	secs	1,5				
Ramp down time	secs	1,5				
Work time	secs	8				
Rest time	secs	20				
Alternating						
Synchronous						
Overall time	599 min					

Programmes

Progr.: 13		Muscle recovery				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	20	20	10		
Mode		MOD	MOD	W/R		
Frequency work	Hz	2-5	5-10	10		
Frequency rest	Hz					
Pulse duration	µS	250-150	250-150	200		
Modulation time	secs	10	10			
Ramp up time	secs			2		
Ramp down time	secs			2		
Work time	secs			10		
Rest time	secs			10		
Alternating						
Synchronous						
Overall time	50 min					

Progr.: 14		Decongestion				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	20				
Mode		WR MOD				
Frequency work	Hz	40				
Frequency rest	Hz	2-5				
Pulse duration	µS	250				
Modulation time	secs	10				
Ramp up time	secs					
Ramp down time	secs					
Work time	secs	5				
Rest time	secs	30				
Alternating						
Synchronous						
Overall time	20 min					

Programmes

Progr.: 15		Pain therapy				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	15	15			
Mode		BST	Cont			
Frequency work	Hz	2x9 (80 Hz)/sec	110			
Frequency rest	Hz					
Pulse duration	μ S	150	175			
Modulation time	secs					
Ramp up time	secs					
Ramp down time	secs					
Work time	secs					
Rest time	secs					
Alternating						
Synchronous						
Overall time	30 min					

- BURST = burst-mode
CONT = continuous mode
MF = modulating mode (Frequency and/or pulse range change automatically in accordance with a pre-given rhythm.
W/R = intermittent operating and resting mode

Modulation time = exponential

For example: Modulating time 10 seconds at a pulse range of 250 – 150 μ s means:

Start at 250 μ s – exponential decrease (first fast then slow) down to 150 μ s in five seconds, then increase (first slow then fast) back to 250 μ s. This complete cycle lasts for 10 seconds.

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To achieve the best possible effect, one of the electrodes should be put on the motor point of the respective muscle if possible. Try to find out this point by slightly moving the electrode around. The other electrode is fastened on the remote muscle end. However, you also obtain a good effect if you put the electrodes respectively on the opposite muscle ends.

While stimulating, please keep this device in reach, to be able to finish the treatment at once, if current becomes uncomfortable.

Please follow in every case absolutely the references to page 32-33 of this manual.

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In most cases it is very easy to find the suitable points for the placing of the electrodes. It is frequently meaningful to take the hurting place to the middle between the two electrodes. We nevertheless recommend to coordinate the points of electrode-placement with the doctor or physiotherapists if possible. The distance between the electrodes shouldn't be much more then 20 cm.

You see some examples of possible placement-points on the following pages.

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Copyright 2004